

Knee Support

Experience unparalleled relief and recovery with our Knee Support, specifically designed for individuals dealing with ACL, PCL, or meniscus injuries, arthritis, tendonitis, and various sports-related knee injuries. A patella stabilizer and side stabilizers ensure optimal knee alignment. The open patella design and breathable material offer superior comfort, making it easier to adhere to your treatment plan as prescribed by your healthcare provider.



Design Highlights:

- Side stabilizers and patella buttress maintain knee alignment aiding in your recovery.
- The open patella & back design alleviates pressure, and provides a more comfortable fit.
- A low-profile, breathable Neoprene construction keeps your knee cool and dry, preventing irritation and discomfort during prolonged use.
- A precision hinge mechanism moves on two independent planes to more accurately mirror the natural movement of the knee.

Indications for Use:

- PCL/ACL Injury
- Arthritis
- Meniscus Tear
- Sprains & Strains
- Knee Instability
- Tendonitis

Knee Support

Fitting Instructions



Application

Step 1. Unfasten all Velcro straps and apply seated with the knee bent. Position the patella pad on the knee cap and align the hinges with the center of the knee joint.

Step 2. Fasten the calf strap first (this strap prevents slipping) and then the top strap. Fasten the popliteal strap. If slipping occurs, check that the calf strap is above the largest part of the calf muscle and tightened securely.

1)



2)



Care Instructions

- Spray with a 50% Alcohol/Water Mixture.
- Rinse and Air Dry.
- Do NOT dry clean, iron, or bleach



get well
urgent care

All products are provided by TLC DME.

For more information regarding product application, functionality, and/or warranties, please call TLC DME Patient Care at 1.855.4TLC DME or email patientcare@tlcdme.com.

All DMEPOS products are covered by CMS, work comp, and commercial insurances. If you have not yet met the deductible associated with your insurance plan, you may have a financial responsibility.