

# TKO FX Splint

This adaptable splint helps you find a comfortable fit and effective support for your hand and fingers when recovering from finger sprains, strains and fractures. The TKO provides the injured finger(s) immobilization needed to comply with the treatment plan designed by your provider and be worn comfortably all day & night.



## Design Highlights:

- One splint acts as a Boxer's Fracture Splint, Radial Gutter Splint, and Middle/Ring Finger Splint.
- Straps are adjustable for fit and comfort.
- Distal buddy strapping eliminates rotation of the fingers.
- Padded stays are malleable and allow for a comfortable and customized fit.

## Indications for Use:

- Boxer's Fracture (ulnar gutter)
- Radial Gutter Fracture
- Middle & Ring Finger Fractures
- Finger Tendon Strains or Sprains, Tendonitis, Arthritis

# TKO FX Splint

## Fitting Instructions



The TKO FX Splint is adaptable to fit your providers treatment plan for a variety of finger sprains and fractures. It is designed to accommodate each of the indications and fitting instructions indicated below. The splint fits either right or left hand.

- Secure the smallest strap first (fingers)
- Secure the middle strap second
- Secure the wrist strap last
- Adjust for proper fit for your indication (below)

### **Boxers Fracture:**

*Left splint*, fits left hand 4th and 5th metacarpals;

*Right splint*, fits right hand 4th and 5th metacarpals.

### **Splinted Middle/Ring Fingers::**

*Left splint*, fits left hand 3rd and 4th metacarpals;

*Right splint*, fits right hand 3rd and 4th metacarpals.

### **Radial Gutter:**

*Left splint*, fits right hand 2nd and 3rd metacarpals;

*Right splint*, fits left hand 2nd and 3rd metacarpals.

### **Care Instructions**

- Hand wash with mild detergent
- Rinse and air dry
- Do NOT use bleach



All products are provided by TLC DME.

For more information regarding product application, functionality, and/or warranties, please call TLC DME Patient Care at 1.855.4TLC DME or email [patientcare@tlcdme.com](mailto:patientcare@tlcdme.com).

All DMEPOS products are covered by CMS, work comp, and commercial insurances. If you have not yet met the deductible associated with your insurance plan, you may have a financial responsibility.